

Int Poster J Dent Oral Med 2001, Vol 3 No 1, Poster 66

# **Elastase Activity in Gingival Crevicular Fluid of Healthy Volunteers**

Language: English

Author(s): Jens Martin Herrmann, José Roberto Gonzáles, Andreas Kleinsteuber, Julia Vonholdt, Jörg Meyle Justus-Liebig-Universität, Giessen, Germany

Date/Event/Venue: 10.03.99-13.03.99

77th General Session of the IADR Vancouver, Canada

### Introduction

Elastolytic enzymes play an important role in the destruction of periodontal tissues. Egelberg (1966) described the influence of gingival and periodontal inflammation on the permeability of the dento-alveolar vessels [4]. Gustaffson related periodontal destruction to the activity of polymorphonuclear neutrophils (PMNs) and their products [5]. The neutrophil acts as the first line of defense against invading pathogens and it is the primary cell present of the acute inflammatory response. Neutrophil dysfunction results in progression of disease and early tooth loss, underlining the primary protective role of PMNs. In addition to elastase, collagenase and cathepsin are other components of the primary or azurophilic granules of neutrophils. Damiano (1988) analyzed the enzyme concentration of azurophilic granula and discovered a high substrate specificity [2]. Great quantities of these enzymes are released during the degranulation process or cell lysis. Janoff (1985), Velvart and Fehr (1987) studied the elastolytic proteolysis of collagen, laminin, fibronectin and proteoglycans [6, 9]. Free elastase is rapidly inactivated by a-1- proteinase-inhibitor. Peroxide radicals usually present in the inflammatory exudation react against this inactivation process [7, 8]. Elastase can therefore also be detected in gingival crevicular fluid of clinically healthy periodontal tissues [3]. Longitudinal studies in patients were performed basically by Cox and Eley (1987), demonstrating an increase in elastasolytic activity in gingival crevicular fluid [1]. Elastase can be recovered from gingival crevicular fluid with Periopapers. Qualitative and quantitative laboratory analysis is possible.

#### Objective

It was the aim of the present study, to establish a micro test system for the detection of elastolytical activity in gingival crevicular fluid. An initial application was tested in experimental gingivitis.

#### **Material and Methods**

#### Calibration of the Periotron 8000®

The Periotron® (Vetter Laborbedarf, Ammerbuch, Germany) was calibrated with double-distilled water and pool serum using a Hamilton® syringe (Hamilton, Darmstadt, Germany) from 0.1 to 0.8  $\mu$ l in volume increments of 0.1  $\mu$ l. The arithmetic mean of 10 measures was calculated and the corresponding calibration was plotted.

#### Development of test and sample buffer

The optimum pH of the assay buffer was found at 8.1. The composition of the buffer used was: Tris-HCl 50 mmol/l (Amersham Life Science, Braunschweig, Germany); sodium chloride 1 mol/l (Merck, Darmstadt, Germany) and Triton-X-100 0.1 % v/v (Sigma, Deisenhofen, Germany). Sample buffer was composed of n-morfolino disodium salt 50 mmol/l (Fluka; Neu-Ulm, Germany); sodium chloride 150 mmol/l (Merck, Darmstadt, Germany); 1.4dithiotreitol 0.1 mmol/l (Merck, Darmstadt, Germany) and Triton-X-100 0.1 % v/v (Sigma, Deisenhofen, Germany) at pH 5.5. (Amersham Life Science, Braunschweig, Germany); sodium chloride 1 mol/l (Merck, Darmstadt, Germany) and Triton-X-100 0.1 % v/v (Sigma, Deisenhofen, Germany); sodium chloride 1 mol/l (Merck, Darmstadt, Germany). Sample buffer was composed of n-morfolino disodium salt 50 mmol/l (Fluka; Neu-Ulm, Germany); sodium chloride 150 mmol/l (Merck, Darmstadt, Germany); 1.4dithiotreitol 0.1 mmol/l (Merck, Darmstadt, Germany) and Triton-X-100 0.1 % v/v (Sigma, Deisenhofen, Germany); at pH 5.5. (Amersham Life Science, Braunschweig, Germany); sodium chloride 1 mol/l (Merck, Darmstadt, Germany); sodium chloride 150 mmol/l (Merck, Darmstadt, Germany). Sample buffer was composed of n-morfolino disodium salt 50 mmol/l (Fluka; Neu-Ulm, Germany); sodium chloride 150 mmol/l (Merck, Darmstadt, Germany); 1.4dithiotreitol 0.1 mmol/l (Merck, Darmstadt, Germany) and Triton-X-100 0.1 % v/v (Sigma, Deisenhofen, Germany) at pH 5.5.

#### **Determination of recovery rates**

A thin space was created between two cover slips and different sample volumes of human leukocyte elastase (Calbiochem, Bad Soden) with known activities (200, 1600 and 3200  $\mu$ U/ $\mu$ I) were introduced and absorbed with a Periopaper®. The volume was determined, then eluted in sample buffer and stored at -22 °C for 48h before enzymatical analysis.

#### **Clinical Study**

A longitudinal, randomized, double blind experimental gingivitis study was performed in 12 healthy male young volunteers, between 22 - 25 years of age. After their written informed consent, a 18 days pre-experimental phase was started with professional prophylaxis and strict control of oral hygiene (days -18, -14, -7, -4). The following clinical parameters were assessed: a Quigley & Hein plaque index (PLI) and the Saxer & Mühlemann papillary bleeding index (PBI). Two contralateral teeth (tT): 13 and 25 were chosen in the upper jaw as test teeth in all volunteers. On day -1 the use of dental floss was stopped. At baseline, all hygiene procedures were completely ceased for the next 18 days. On day 0, 3, 7, 14 and 18, clinical parameters were assessed and GCF was taken from the mesiobuccal site of both tT, placing two Periopapers® for 20 seconds in the margin of the sulcus. Immediately afterwards, the fluid volume was measured in the Periotron 8000® and the two strips corresponding to each tooth were pooled into an Eppendorf® cup with the respective sample buffer.

#### Analysis of GCF samples

After a centrifugation at 5000 rpm for 5 min., 10 µl of sample were added to a microtiterplate containing the assay buffer, together with 50 µl of a 10-3 mol/l solution of the fluorogenic substrate MeO-Succ-Ala-Ala-Pro-Val-7-amino-4-methylcoumarin (Bachem Biochemica, Heidelberg, Germany). The plates were covered with a removable film and incubated for 6 h at 25°C.




Fig. 1: Calibration of the Periotron 8000 a) r<sup>2</sup>double distilled water = 0.999 b r<sup>2</sup>poolserum = 0.998 (polinomial regression 3rd degree) Each sample activity was directly compared to a standard curve of human leukocyte elastase (HLE). Elastolytic activity of the samples was directly calculated using the Biolise® software (SLT-Labinstruments, Crailsheim, Germany). Each sample was analyzed twice.




Fig. 2: Standard Curve of the Elastolytic Activity F.U. = Fluoreszenz units - Activity =  $nU/\mu l$  Excitation: lambda = 390 nm; Emission: lambda = 460 nm

The proteolytic cleavage of the polipeptide substrate liberated AMC. The emitted light (lambda = 390 nm) was quantified at 460 nm in the Fluostar®.




Fig. 3: Sample stability during storage at different temperature

#### Results

Calibration curves of the Periotron 8000® are shown in Fig. 1. An optimum curve fit was obtained with a polynomial regression calculation ( $r^2 = 0.999$ ). The fluorogenic substrate MeO-Succ-Ala-Ala-Pro-Val-7-amino-4-metylcoumarin can be used for the diagnosis of elastolytic activity. A standard curve of human leukocyte elastase is shown in Fig. 2.




Fig. 4: Mean assay recovery rates (n=50)

The sample buffer allows a storage time up to 48 h (at a pH of 5.5). Better stability of elastolytic activity was observed when samples were kept at -80 °C or at -22 °C (Fig. 3). It is interesting to note that stored at +4 °C, the activity began to decrease after 12 h. After 72 h they presented only 84 % of their initial level.

The recovery rates of this method were about 95 % for all three activities. Activity choice was done in relation to expected clinical activity ranges (Fig. 4).




Fig. 5: PLI = plaque index and PBI = papilary bleeding index (during the experimental phase of the study)

Test limit was about 1  $\mu$ U/ $\mu$ l. However, samples could be recovered from a specific site of the tooth and could be analyzed even in small concentrations of 0.1  $\mu$ l.

Elastolytical activity was analyzed in all volunteers, even without presenting gingivitis. Fig. 5 shows the median PLI and PBI. At initial examination, both clinical parameters exhibited similarly low scores.

After baseline, PLI increased continuously until day 18 (about 2.8 fold). An increase of the PBI was observed from day 7 to day 18 (about 1.0). Median elastolytic activities of the test teeth are shown in Fig. 6. The median activity in gingival crevicular fluid was 452  $\mu$ U/ $\mu$ l at baseline, 708  $\mu$ U/ $\mu$ l at day 3, 1269  $\mu$ U/ $\mu$ l at day 7, 1325  $\mu$ U/ $\mu$ l at day 14 and 1507  $\mu$ U/ $\mu$ l at day 18.




Fig. 6: Changes of elastolytic activity during experimental gingivitis in 12 human volunteers (red = two dropouts)

#### **Discussion and Conclusions**

Much effort has been directed since the past decade in the identification of host and bacterial enzymes in crevicular fluid.

Detection of minimal elastolytical activity was possible, even considering that the alkaline assay pH cannot be compared to in vivo pH as found in inflammation. Our assay reflects levels of the total elastolytic activity produced by PMN lysis and dissolved total elastase. The samples were maintained under conditions appropriate for minimizing loss of activity. Our study shows that elastase activity in GCF increases strongly with the development of experimental gingivitis. Although several different assay procedures to measure elastolytic activity have been published, only in a very few of them, the quality of the assay procedure was checked. In our system recovery rates were above 95% and were considered acceptable for use of the assay. The elastolytic activity, that was found in GCF of volunteers, could be due to the different assay procedure.

Further advantages of this method are the possibility to analyze specific sites at single teeth, optimal utilization of the reagents, as well as analysis of up to 80 samples per reading.

### Bibliography

- 1. Cox, S. W., Eley, B. M.: Preliminary studies on cysteine and serine proteinase activities in inflamed human gingiva using different 7-amino-4trifluoromethyl coumarin substrates and protease inhibitors. Arch Oral Biol, 1987, 32, S. 599 - 605
- Damiano, V. V., Kucich, U., Murer, E., Laudenslager, N., Weinbaum, G.: Ultrastructural quantitation of peroxidase- and elastase-containing granules in human neutrophils. Am J Pathol, 1988, 131, S. 235 - 245
- Darany, D. G., Beck, F. M., Walters, J. D.: The relationship of gingival fluid leukocyte elastase activity to gingival fluid flow rate. J Periodontol.1992, 63, S. 743 747
- 4. Egelberg, J.: The blood vessels of the dento-gingival junction. J Periodontal Res, 1966, 1, S. 163 179
- Gustafsson, A., Asman, B., Bergstrom, K.: Elastase and lactoferrin in gingival crevicular fluid: possible indicators of a granulocyte-associated specific host response. J Periodontal Res, 1994, 29, S. 276 - 282
- 6. Janoff, A.: Elastase in tissue injury. Annu Rev Med, 1985, 36, S. 207 216
- 7. Kowashi, Y., Cimasoni, G., Matter, J.: Collagen breakdown by gingival collagenase and elastase. Experientia, 1980, 36, S. 395 396
- 8. Ohlsson, K.: Alpha1-antitrypsin and alpha2-macroglobulin. Interactions with human neutrophil collagenase and elastase. Ann N Y Acad Sci, 1975, 256, S. 409 421
- 9. Velvart, M., Fehr, K.: Degradation in vivo of articular cartilage in rheumatoid arthritis and juvenile chronic arthritis by cathepsin G and elastase from polymorphonuclear leukocytes. Rheumatol Int, 1987, 7, S. 195 202

This Poster was submitted on 08.02.01 by Jens Martin Herrmann.

#### **Correspondence address:**

Jens Martin Herrmann University of Giessen Dept. of Periodontology Schlangenzahl 14 35392 Giessen Fax: +49 (0)641 / 9946189 # 3242

## Elastase Activity in Gingival Crevicular Fluid of Healthy Volunteers J. M. HERRMANN<sup>\*</sup>, J. R. GONZALES, A. KLEINSTEUBER, J. VONHOLDT, J. MEYLE

Alter is constitution at 200 week milled to is microther

Justus - Liebig - Universitaet, Glessen, (Germany)

and at 8,1 The HCT SD menal A emory: sodium

#### ARSTRACT

The role of coll agen the Not only the

NSO INC

OBJECTIVE

MATERIAL AND METHODS

It was the win of the present dudy, to establish a mice test and we for the detection of eliadorytical advicts in graphic previouslie full, an initial application was to deal in experimental graphics.

Calibration of the Pendotro ISSS  $\vec{0}$ The Pendotro I (Vetter Laboritador / Ammado Jah, Germany) and their deviation of Vetter Laboritador / Ammado Jah, Germany and end of "penge Ham Itoo, Dann ded, Germany Teori (Tho Gal Jah Seduna Transmontost OI (Jah, The attras decisionari measure was and adulted and the corresponding calibration sample data

an descripted period of this were end under any applicated electra as in the GCP of all w electrolytics octivity was 652 µUUµ. In or or ignoups this activity can be explained set if in this tend are note to ble on the espectro-center of the ender on the to ble on the espectro-

stepflight fill placing two he mangin of the suicus, volume was measured in For for strength



RESULTS C mitmation curves of the Periotion 3000' are shown in Fig. 1. An optimum curve III was obtained with a polynomial regie score Ale-Ma-Pic-Val-7-emitte le diagnosis of desitant



CONCLUSIONS

 T. & J. W.M. May, A.M. Paterson, J. Isasam ginglast using although 7-mails and toward-on Tax 32, 400 (1997) tari comilei n in 10, 20, 10, 100 Cased y THE OTHER OF a Revenue Aur 1, 1000 times BITTY APPENDICATION CONTRA ranansial artein and tarein ine transfer Alexandri in 7. March 11

containing the account of the fluo



2

cleavage of the p that light (% = 3.50



in Fig. 2. The sample built (at a pH of 55), observed when sam



cletical add vityrainges (Fig.4.). Ted. Will was about 1 µli (µl. H diffie even in as all concentrations of 0.1 µl. Elattorytical activity was maligzed in all volunt eers, even without Eletator(#doi) activ Evves and stabil in elivatives; even valicut area may grapp de . Fgl. 5 pay offen elevative random elivative for the elevative and the balance area para transmission and an elivative provember of the second elivative and the elivative (a 10). Allo introduce d'hiele Bi valeo losovice el en de y? To lety 1 (a 10).

